



Straight and Kurved amniotic sac membrane perforator

A Straight amniotic membrane perforator tool with sharpened point

- High visibility purple colour
- Sterile (ETO sterilised)
- Single use
- Disposable

A Kurved amniotic membrane perforator tool with sharpened point

The Kurve Hook perforator is gently curved to follow the natural contours of the birth canal.

Features

- High visibility purple colour
- Sterile (ETO sterilised)
- Single use
- Disposable

Technical Specifications

- Approximately 28cm long
- Polypropylene
- Extra sharp point
- CE certified

Ordering Information

Straight

Part Number	SH-1
Quantity	100 units per case
Description	Straight Hook
Delivery	Usually ex-stock

Ordering Information

Kurved

Part Number	KH-1
Quantity	100 units per case
Description	Kurve Hook
Delivery	Usually ex-stock