



Spackman Cannula

The unique device with an atraumatic elastomeric tip and an adjustable cervical cone that ensures an effective seal during use.

The Spackman Cannula is unique in having both an adjustable cervical cone and an atraumatic elastomeric tip.

The adjustable cervical cone allows the Cannula to be 'made to fit' the depth of the uterus which ensures an effective seal and reduces the risk of dye leakage. The atraumatic elastomeric tip is designed to reduce the risk of perforation to the uterus.

The Spackman Cannula can be purchased separately or as part of the DiOne procedure pack.

Ordering Information

| | |
|-------------|-------------------|
| Part Number | SPK-01 |
| Quantity | 25 units per case |
| Description | Spackman Cannula |
| Delivery | Usually ex-stock |