



Cyclopes

Self contained, self illuminating proctoscope

A unique proctoscope with its very own internal light source designed to provide a shadow-free environment during rectal examinations.

Cyclopes is a single patient use, self illuminating proctoscope. The white body and innovative internal light source provides a shadow-free environment ensuring that there is no need for any additional light source or power for operation.



The unique design of Cyclopes is an innovative light source placed at the base of the obturator. Removing the obturator automatically switches the light source on, replacing the obturator automatically switches the light off. When you have finished using the proctoscope the light source can be removed. The circuit board and battery can be disposed of independently based on local hospital/surgery protocol.

The tapered front of the obturator allows for gentle insertion into the anal canal. The handle is ergonomically designed to fit comfortably into your hand and provides full support in the manipulation of the proctoscope into position.

Cyclopes comes sterile packed and ready for immediate use.

Ordering Information

Part Number	CLP-01
Quantity	25 units per case
Description	Self illuminating proctoscope
Delivery	Usually ex-works

Technical Specifications

Moulded from medical grade ABS, white	
Tapered insertion end length 20.75mm, taper 24.2°	
Insertion length 50mm, width 25mm	
Visualisation diameter 20.2mm, area 320.31mm	
Handle length 92mm, width 26mm, depth 10mm	
Total length of Cyclopes	145.7mm

Electrical Specifications

Light source: white LED 5mm diameter, 1200 med @ 15°	
EMC: EN 60601-1:1990	
MES to EN 60601-1:1990	
LED intensity to EN 60825	
Manufactured and tested to ISO 13485:2003	