



# The Complete Fetal Scalp Blood Sampling procedure pack

Including the LiNA™ Fetal Blood Sampling Wand

OBSTETRICS

The Complete Fetal Scalp Blood Sampling (FSBS) procedure pack contains everything necessary to complete a fetal blood sample, without the use of petroleum jelly.

The FSBS procedure pack contains a LiNA™ single handed Fetal Blood Sampling Wand, four heparinised capillary tubes, capillary caps, mixing wires and universal blood gas analyser adapters. The adapters will fit the majority of commonly available blood gas analysers. The tray is wrapped in a sterile surgical drape for easy dispensing.

The FSBS procedure pack contains three extra capillary tubes if more samples are required.

The procedure pack can also be purchased with the ELA self-contained, self-illuminating amnioscope.

#### Capillary Tube Specifications

The capillary tube is sodium heparinised 240IU/ml

Inner diameter	1.4mm
Outer diameter	2.0mm
Length of capillary tube	85.0mm
Capillary volume	minimum 130µl

#### Principals of Operation

- 1 The knife at the end of the sampler is pushed out to penetrate the scalp and then retracted for safety.
- 2 The capillary tube is used to collect a blood sample.

#### Ordering Information

Part Number	4201-AC
Quantity	5 units per case
Description	Complete Fetal Scalp Sampling Kit
Delivery	Usually ex-stock

#### Technical Specifications 4201-AC

The tray contains the following\*

1	LiNA™ FBS 4201 Fetal Blood Sampler
6	Long cotton tipped buds
4	Heparinised capillary tubes
8	Capillary caps
4	Mixing wires
4	Universal blood gas analyser adapters with clot buster
1	Sterile surgical drape
A resting space for ELA-01 amnioscope	

\* Products are subject to change

Petroleum jelly free fetal blood sampling

[www.bmmedical.co.uk](http://www.bmmedical.co.uk)

BridgeMaster Medical Limited, Unit 4, Ashville Trading Estate, Suffolk Way, Abingdon, Oxfordshire OX14 5JX, United Kingdom  
Tel +44 (0)1235 858290, Fax +44 (0)1235 530033, Email [info@bmmedical.co.uk](mailto:info@bmmedical.co.uk)